Snelgids

© Sig vzw
versie 22 oktober 2021

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patient sessions
your patient sessions

3. Neem de module SKO.DDK af

experiments
record patient session

4. Haal resultaten op

patient sessions
progress & results

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Add a new patient with the minimum number of details. You can add more details later under 'your patient'.

Patient profile

surname (*): Peeters
first name (*): Jan
date of birth: year = 2016 month = 04 day = 02

1

2

Patient "Peeters Jan" is successfully added.

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patient sessions **your patient sessions**

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

1 Select a template or session:

your patient sessions Select a template or session: SKObasis, SIG, vzw, ga

Session properties date of creation (yyyy-mm-dd): 2021-05-09
description (maximum 32 characters): vooorbeld
selected patient: SKObasis, SIG, vzw, ga

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

3 Select a template or session: 2021-09-30, Peeters, Jan

your session templates Select a template or session: diadochokinese

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

2 Select a template or session: SKO, ddk1[SIG,vzw,ga]

your patient sessions Select a template or session: SKO, ddk1[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

your session templates Select a template or session: SKO, ddk2[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

your session templates Select a template or session: SKO, ddk2[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

your session templates Select a template or session: SKO, ddk2[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

your session templates Select a template or session: SKO, ddk2[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

your session templates Select a template or session: SKO, ddk2[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

Supported actions Add new patient sessions and update or delete formerly defined patient sessions

your session templates Select a template or session: SKO, ddk2[SIG,vzw,ga]

Session properties date of creation (yyyy-mm-dd): 2021-09-30
description (maximum 32 characters): diadochokinese
selected patient: Peeters, Jan

The patient session was successfully added.

SKO – © SIG vzw – www.SKOTest.be

3. Neem de module DDK af

experiments **record patient session**

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test microphone try exercise record patient session

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Supported actions test microphone try exercise record patient session

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Cyclus 2 & 3

4. Haal resultaten op

patient sessions **progress & results**

Supported actions Survey of your active patients. Per patient, you see the number of sessions t completed by the patient. Click on a **sessions** button to consult more details for the corresponding pati

patient name	nr of sessions defined / completed	sessions
But Jeanne	1 / 1	sessions
But Jerom	1 / 1	sessions
Caro Lientje	3 / 3	sessions
DDK Joris	7 / 7	sessions
john F	0 / 0	sessions
leon jef	1 / 1	sessions
M29190 VVA	1 / 1	sessions
More Thomas	1 / 1	sessions
Name Louise	1 / 0	sessions
Peeters jan	1 / 1	sessions
Pek Pietje	3 / 3	sessions
Slagers Georges-Louis	2 / 2	sessions

Supported actions Survey of the sessions that were created for patient Peeters jan. Click on a **more** button (if any) to consult more details about the self-evaluations, the automatic a particular session.

one level back	patient session + description	session cycles (due / done)	session analyses (due / done)	session ratings (open / busy / done)
your patient sessions	2021-09-27 (diadochokinese)	3 / 3	6 / 6 more	0 / 0 / 0

Automatic analysis outputs obtained for a particular session. To print the report: put the cursor in this window, right click on the mouse and select the PRINT opt

one level back

Peeters jan (session: 2021-09-27 diadochokinese)

AMR= papapapapapa..... (14 syllables)
SMR= potteke potteke.... (15 syllables)

Supported actions Automatic analysis of speech recordings

cycle	analysis method	select figures	analysis result	comment
3	AMR.PRAAT (spreeksnelheid)		SR = 6.12 syllables/sec (14 syllables and 0 pause intervals)	
	SMR.PRAAT (spreeksnelheid)		SR = 8.52 syllables/sec (14 syllables and 0 pause intervals)	
2	AMR.PRAAT (spreeksnelheid)		SR = 4.19 syllables/sec (13 syllables and 0 pause intervals)	
	SMR.PRAAT (spreeksnelheid)		SR = 6.15 syllables/sec (16 syllables and 0 pause intervals)	
1	AMR.PRAAT (spreeksnelheid)		SR = 4.57 syllables/sec (14 syllables and 0 pause intervals)	
	SMR.PRAAT (spreeksnelheid)		SR = 6.09 syllables/sec (15 syllables and 0 pause intervals)	

SKO - © Sig vzw - www.SKOTest.be